

Laser Choices – Overview

	Picasso PowerDesktop	QuickLase PowerLase Cordless
	7w	4w
Battery / Mains	Mains	Battery Long life / fast charge
Colour Touch Screen Multi-languages	✓	No need
Fibre optic housing	✓	No need
Fibre type	Strippable/Disposable tips	Disposable tips
Fibre Optic Size	200-400 micron	200-400 micron
Video Procedures	✓	✓
810nm wavelength	✓	✓
Ultimate cutting	✓	✓
Ultimate Coagulation	✓	✓
Easy to use	✓	✓
Continuous Wave mode	0.1w – 7w	auto power
Pulsed Wave mode	✓	auto set
Sterilisable hand pieces	✓	chemi wipe
Foot pedal operated	✓ cordless/corded	Index finger touch
Customisable pre-sets	8	2 sterilisation/cutting

Perio, Endo, Ortho	✓	✓	✓
Hygiene	✓	✓	✓
Spot Whitening	✓	✓	✓
Troughing	✓	✓	✓
Curettage	✓	✓	✓
Frenectomy	✓	✓	✓
Ulcers	✓	✓	✓
Crown Lengthing	✓	✓	✓
Gingivectomy	✓	✓	✓
Papillectomy	✓	✓	✓
Implant Recovery	✓	✓	✓
Sulcular debridement	✓	✓	✓
Soft tissue management	✓	✓	✓
Less traumatic	✓	✓	✓
Less Anaesthesia	✓	✓	✓
Fast Recovery	✓	✓	✓
Control bleeding	✓	✓	✓
Sterilisation	✓	✓	✓
Reduced fatigue & Stress	✓	✓	✓
Less training	✓	✓	✓
Decreased practice time	✓	✓	✓
No collateral damage	✓	✓	✓
Tissue Welding (No suture)	✓	✓	✓
Increase in productivity	✓	✓	✓