


QuickLase™

Laser specialist since 1992



3w (810nm)
5w (810nm)
8w Dual (4w 810nm & 4w 980nm)


W. 15cm (6") L. 23cm (9") H. 10cm (4") Weight 1.6kg (3lb 8oz)

Many dentists are discovering the benefits of using a laser for routine treatments. The need for sutures is reduced with healing time and bleeding dramatically reduced.

A further benefit is that many treatments can be carried out in one visit - convenient and quick.

A choice of soft tissue diode lasers with different wavelengths; 810nm and 810nm+980nm for ultimate cutting and coagulation.

Ultimate Cutting with Ultimate Coagulation



Research has shown that the most efficient soft tissue laser utilises two wavelengths : 810nm and 980nm. They work in unison to optimise coagulation and ablation.

- Mains and battery power options
- 2 Year warranty
- Manufactured and serviced by Quicklase in the U.K.
- Buy direct and save thousands



10w (5w 810nm + 5w 980nm)
20w (10w 810nm + 10w 980nm)
Autoclaveable FibreCaddy™

W. 16cm (6 1/4")
L. 26cm (10 1/4")
H. 12cm (4 3/4")
Weight 1.8kg (4lb)

QuickLase™

www.quicklase.com - info@quicklase.com

+44 (0) 1227 780009

FDA CE ISO
Certified dental procedures

