


Laser Accessories

Fibre tips handpiece



Handpiece for fibre tips, therapy and biostimulation as well as whitening tips.

Product Code: QWDISHAND

TMJ Tips




Biostimulation Tips



Whitening Tips



Bendable fibre tips



For quick and easy use with the bendable fibre tips handpiece, for gum cutting, endo, perio and therapy biostimulation. Comes in 400 or 200 microns colour-coded, can be used as disposable or autoclavable you can specify on order notes.

Product Code: QWDISTIPS



Tip	Diameter[μm]	Length[mm]
5mm	300	5
5mm	400	5
10mm	300	10
10mm	400	10
10mm	600	10
15mm	200	15



Uninterruptible Power Supply



Advanced Power Protection

Your laser will be protected on a stable power network. This UPS stabilizes and refines the power supply, protecting the laser from power disturbances.

Product Code: QWUPS



Functional Ease-of-Use

Front panel LED indicators with audible alarms to easily tell the functional status.

AVR

Automatic Voltage Regulation

Instantly corrects incoming utility power without utilizing the battery, saving the battery for when it is needed most.