


Quicklase Laser DentaLase model Choices – Overview

	QuickLase 5w 	QuickLase 6w or 12w Dual 	QuickLase 12w Dual plus 	QuickLase 6w 1064 
Wavelength	810nm	810nm + 980nm at the same time or individual and adjustable	810nm + 980nm at the same time or individual and adjustable	1064nm
Clinical procedure pre-sets	20 per-set procedures	20 pre-set procedures	24 pre-set procedures with 5 short cuts	16 pre-sets
Colour coded Touch Screen Multi-languages	✓ 3" screen 6 languages	✓ full colour, 4" screen 6 languages inc Arabic and clinical pictorials	✓ full colour, 6" screen 6 languages inc Arabic and clinical pictorials	✓ 3" screen 6 languages
Fibre optic housing	✓	✓	✓ With coloured LED	✓
Strippable fibre as standard with disposable tips - more economical <i>(non-strippable fixed fibre with bendable disposable tips - available on request)</i>	✓	✓	✓	✓
Fibre Optic Size	200-400 micron	200-400 micron	200-400 micron	200-400 micron
Bleaching Procedure Bleaching prism included	✓	✓	✓	✓

Therapy Procedure Therapy prism included	✓	✓	✓	✓
Ultimate coagulation & lasing/cutting Dual technology		✓	✓	
Easy to use	✓	✓	5 clinical shortcuts	✓
Continuous Wave mode	0.1w – 5w	0.1w – 12w	0.1w – 12w	0.1w – 6w
Pulsed Wave mode 20,000 Hz	20, 30, 50Hz & adjustable	20, 30, 50Hz & adjustable 20,000 Hz	20, 30, 50Hz & adjustable 20,000 Hz	20, 30, 50Hz & adjustable
Sterilisable All-in-1 handpiece	✓	✓	✓	✓
Complete with accessories Fibre, stripper, cleaver, hand-piece, caution yellow sticker, 25 disposable tips, USB stick training and docs Whitening and therapy prisms	✓	✓	✓	✓
Wired Foot pedal operated (Wireless available on request)	✓	✓	✓	✓
3 protective glasses	✓	✓	✓	✓
Comes with Bleaching & Therapy prisms	✓	✓	✓	✓
Comes with patients leaflets & posters	✓	✓	✓	✓
Training	On-line & USB stick	On-line & USB stick	On-line & USB stick	On-line & USB stick
BILD hands on training	available	available	available	available

QuickLase™ DUAL



Some of the Clinical procedures

	5w 810nm	6w Dual 810+980nm	12w Dual 810+980nm	12w Dual + d810+980nm	6w 1064nm
Perio, Endo, Ortho LLLT	✓	✓	✓	✓	✓
Hygiene	✓	✓	✓	✓	✓
Spot Whitening	✓	✓	✓	✓	✓
Troughing	✓	✓	✓	✓	✓
Curettage	✓	✓	✓	✓	✓
Frenotomy	✓	✓	✓	✓	✓
Ulcers	✓	✓	✓	✓	✓
Crown Lengthening	✓	✓	✓	✓	✓
Gingivectomy & de-pigmentation	✓	✓	✓	✓	✓
Papillectomy	✓	✓	✓	✓	✓
Implant Recovery	✓	✓	✓	✓	✓
Sulcular debridement	✓	✓	✓	✓	✓
Soft tissue management	✓	✓	✓	✓	✓
Less traumatic	✓	✓	✓	✓	✓
Less Anaesthesia	✓	✓	✓	✓	✓
Fast Recovery	✓	✓	✓	✓	✓
Control bleeding	✓	✓	✓	✓	✓

Sterilisation	✓	✓	✓	✓	✓
Reduced fatigue & Stress	✓	✓	✓	✓	✓
Less training	✓	✓	✓	✓	✓
Decreased practice time	✓	✓	✓	✓	✓
No collateral damage	✓	✓	✓	✓	✓
Tissue Welding (No suture)	✓	✓	✓	✓	✓
Increase in productivity	✓	✓	✓	✓	✓
Better Endo sterilisation	✓	✓	✓	✓	✓