

QUICKLASE LASERS



5w
810nm

QWLASER5
(3.5" Screen)



6w/12w Dual
810 + 980nm


QWLASER6/12D
(4" Screen)



12w Dual Plus
810 + 980nm

QWLASER12DP
(6" Screen)

- 7 Key applications: Surgery - Implants - Coagulation - Perio - Endo - LLLT Therapy - Whitening
- Choice of 810nm, 810nm+980nm for better coagulation and ablation
- Presets for most used clinical procedures
- Built-in fibre caddy for fibre protection, fibre 400 & 200 micron option
- 3.5" Colour coded touch screen with 20 Presets (5w)
- 4" Full Colour coded touch screen with 20 Presets (6w/12w Dual)
- 6" Full Colour coded touch screen with 24 presets and 5 Shortcuts (12w Dual Plus)
- 6 languages including Arabic
- Adjustable aiming beam with instant on/off function
- Built in auto step lase
- 3 Pulse modes adjustable & Continuous mode cutting
- Preset 30 second timer for whitening & therapy
- Overheat protection with automatic reliable cooling system
- 100% EMC protective case
- Complete with accessory package and online training
- Laser Hands-on courses also available with BILD Academy
- See our website for all additional accessories





No expensive
tips



No battery
No power loss



Pain free
laser dentistry



QuickLase™
+44 1227 760649
www.quicklase.com

